

Soft Ionization Spectrometry (SIS):

***One Step Closer to a “Magic Box” for
Complete Bev-Gas Testing – AND Beyond***

Presented at the ISBT Conference – San Antonio - 2007



Airborne Labs International

Don Pachuta, Ph.D. & Joe Angeloni

Background - Review

Industrial Sources of CO₂

+ Unique Distribution Chain

- Natural Wells
- Fuel Combustion
- Bio-Fermentation
- Ammonia Production
- Chem-Process Off Gas



- Feedgas Source – Producer – Transport – Regional Depots – Transport – Storage – Production

Many Potential Impurities

ISBT Bev-Gas Quality Challenges

- As beverage production grows, *continuous* quality testing becomes v important
- CO₂ Quality Assurance is a complex technical challenge
- All dream of a “Magic Box” that can quickly monitor the *entire ISBT Impurity List*

ISBT CO₂ Purity Guidelines

"The List"

| | | |
|--|---|----------------------|
| CO₂ Purity | 99.9% min | Process |
| Moisture (H₂O) | 20 ppm v/v max | Process |
| Oxygen (O₂) | 30 ppm v/v max | Sensory |
| Carbon Monoxide (CO) | 10 ppm v/v max | Process / Regulatory |
| Ammonia (NH₃) | 2.5 ppm v/v max | Process |
| NO / NO₂ | 2.5 ppm v/v (ea) max | Regulatory |
| Non-Volatile Residue (NVR) | 10 ppm x/x max | Sensory |
| Non-Volatile Organic Residue (NVOR) | 5 ppm w/w max | Sensory |
| Phosphine (PH₃) | 0.3 ppm v/v max | Regulatory |
| Total Volatile Hydrocarbons (THC) & TNMHC | 50 ppm v/v max – including 20 ppm v/v TNMHC | Sensory |
| Acetaldehyde (AA) | 0.2 ppm v/v max | Sensory |
| Aromatic Hydrocarbon Content (AHC) | 20 ppb v/v max (as benzene) | Regulatory |
| Total Sulfur Content (TSC*) less SO₂ | 0.1 ppm v/v max | Sensory |
| Sulfur Dioxide (SO₂) | 1 ppm v/v max | Sensory |
| Odor of Snow | No foreign odor | Sensory |
| Appearance in Water | No color or turbidity | Sensory |
| Odor & Taste in Water | No foreign odor/taste | Sensory |
| Supplemental HCN / Vinyl Chloride | None detectable–best method | Report |

Where Does SIS Technology *Fit-In* vs Traditional ISBT CO₂ Analyzers?

- CO₂ Purity: Zahm/Nagel or GC/MS

- Moisture: Dew Cup or Hygrometer

- CO DT, FID, NDIR

- NH₃ DT, GC/MS, NDIR

- Phosphine Wet Chem, DT

- THC/TNMHC THA / GC-FID



- Acetaldehyde DT, GC PID, FID

- Aromatic HC BTEX

DT, GC PID, NDIR

- HCN/VCL DT or Wet Chem


- SO₂ or TRS DT, SCD, CLD



What Is Soft Ionization Spectrometry?


- SIS is a Mass Selective Analyzer:
- CO₂ Sample is **continuously** injected into a **Low Energy Ion Beam** of Hg⁺ or Xe⁺ or Kr⁺ = soft ionization gas (SIG)
- Sample impurities (analytes) softly fracture into **only a few** characteristic ion fragments
 - Fragments are **size-separated** by a scanning RF field. Software-selected analyte ions are counted by a pulse detector
 - SIS Display = “**impurity hit list name + ppm level**” — total analysis time = a few seconds
- **So, Tell me more about how this SIS thing works**

The “Soft Ionization” Process




The "Soft Ionization" Process


In Slow Motion



SIS Analyzer System



SIS Analyzer – Under-the-Hood




SIS - a “New Application” of an Established Technology

- SIS is a relatively “mature” technology
- Established Applications:
 - Environmental Air Testing
 - Engine Exhaust / Catalyst / Diesel Studies
 - Medical e.g. Disease Diagnosis by Patient Breath Analysis
- European + Emerging **USA/Canada** CO₂ Supplier QC
 - **PET bottle** recycling & mfg quality testing*
 - Gaseous Headspace Food & Beverage Testing (**Fragrance ID*** & Brand Profiling*)
 - Other Applications where **many** analytes might exist for only milliseconds & must be **quickly** ID'd & quantified


* **AND Beyond???**

PET Bottle Recycling-High Speed QC Testing*

- 1 Machine protection without rejection at the infeed (machine stop with horizontal containers)
- 2 Height detection with laser PE sensor assembly (optional)
- 3 Cap recognition (optional)
- 4 Label monitoring (optional)
- 5 Rejection at machine infeed (optional)
- 6 Mass spectrometer
- 7 Colour spectrometer
- 8 Rejection of contaminated containers
- 9 Rejection monitoring of test containers
- 10 Collecting container




- 1 Neutral air is blown into the container.
- 2 The mixture of air and gas is drawn off and examined for contaminations.




Non-contact sampling

***AND Beyond – for another day**

Fragrance Mass Moment ID Profile*



= Apple juice, peach juice, pear juice, grape juice, apricot juice, tomato juice, sour cherry juice, currant juice, water



*AND Beyond – for another day


So, How might SIS work in the “*Real World*” of Bev-Gas Testing?

Our R&D Objectives & Current Focus

- Develop SIS Methods (SIG, AM, Scan, Count) for key ISBT-listed CO₂ impurities
- Evaluate: LDR / Sensitivity / Precision / Accuracy & DL vs ISBT Method Guidelines
- *Identify potential interferences* + their influence + develop corrective actions
- Evaluate: Analysis - Calibration Strategy + SIS Ops & Maintenance Requirements



SIG & AM Selection Process

Sometimes its "easy"



SIG & AM Selection Process – The “Frag Factor”

(Sometimes its not-so-easy)




Basic SIG & AM Studies

Needed to ID potential Interferences for design of

Appropriate Software & Calibration Routine

Corrective Measures



SIS - Aromatic Hydrocarbon (AHC) Data

(ISBT CO₂ Limit = 20 ppb v/v as Benzene + report all other target AHC's)


Our Results

| AHC | SIG | AM | LDR ppb v/v | DL ppb v/v | Interference (s) |
|--|-----------|------------|-----------------|---------------|--|
| Benzene | Hg | 78 | LT 0 - 1000+ | LT 1 | None found-to-date |
| Toluene | Hg | 92 | LT 0 - 1000+ | LT 1 | None found-to-date |
| Ethyl Benzene + o,m,p Xylenes | Hg | 106 | 1 -1000+ | LT 5 | None found-to-date Note: EBz + Xylenes measured as Total |


Speciation Information B-T-EX

Use of Benzene Cal Gas for Indirect Cal of TEX works v. well

SIS Data: B-T-EX in CO₂ @ 28 ppb v/v level



SIS Data: B-T-EX in CO₂ @ 15 ppb v/v level




SIS - Volatile Sulfur Compounds (VSC) Data

ISBT Limit = 0.1 ppm v/v as TSC* = sum of VSC target impurities excluding SO₂, SO₂ = 1 ppm v/v max


Useful Speciation information achieved

| VSC | SIG | AM | LDR ppb v/v | DL ppb v/v | Interference (s) |
|------------------|-----|----|----------------|---------------|---|
| H ₂ S | Hg | 34 | 0 – 1000+ | LT 5 | PH ₃ @ AM = 34 (Note: PH ₃ = Rare) |
| COS | Xe | 60 | 0 – 1000+ | LT 5 | None found-to-date |
| SO ₂ | Xe | 64 | 0 – 1000+ | 5 | None found-to-date |
| MeSH | Hg | 48 | 0 – 1000+ | LT 5 | None found-to-date |
| CS ₂ | Hg | 76 | 0 – 1000+ | LT 5 | None found-to-date |
| DMS + EtSH | Hg | 62 | 0 – 1000+ | LT 5 | Vinyl Chloride (VCL) AM = 62, correct using AM = 64 isotope (Note: VCL = rare) Note: Measures Total DMS + EtSH |


SIS Data: VSC in CO₂ @ 110 ppb v/v level



SIS Data: VSC in CO₂ @ 60 ppb v/v level



SIS Data: VSC in CO₂ @ 13 ppb v/v level



SIS - Volatile Oxygenated Compounds (VOX) Data


ISBT Limit = 0.2 ppm v/v as Acetaldehyde (AA) + Speciation information needed for THC/TNMHC calcs etc.

| VOX | SIG | AM FR = Frag | LDR ppm v/v | DL ppm v/v | Interference (s) |
|-----------|-----|------------------|----------------|---------------|--|
| AA | Hg | 43 ^{fr} | 0 - 5+ | LT 0.05 | Acetone = Major INT, i-PA = INT, i-Butanes = INT, n-Butane = minor INT |
| DME | Hg | 46 | 0 - 5+ | LT 0.05 | No effect on AA, NO ₂ = INT |
| MeOH | Xe | 31 ^{fr} | 0 - 5+ | LT 0.05 | No effect on AA, EtOH / nPA = INT |
| EtOH | Hg | 45 ^{fr} | 0 - 5+ | LT 0.05 | No effect on AA, DME = minor INT |
| N-Butanol | Hg | 56 ^{fr} | 0 - 5+ | LT 0.05 | No effect on AA |
| Acetone | Hg | 58 | 0 - 5+ | LT 0.05 | Major AA effect, Butanes = INT |
| EtOAc | Hg | 70 ^{fr} | 0 - 5+ | LT 0.05 | No effect on AA |


Other VOX impurities studied to date: n / i -Propanols, i -Butanols, Ethylene Oxide

SIS Data: VOX in CO₂ @ 2.5 ppm v/v Acetone / 0.00 AA – INF Corrected


Key to Success is having accurate Interference Correction Data




SIS Data: VOX in CO₂ @ 2.5 ppm v/v Acetone + AA – INF Corrected




SIS Data: VOX in CO₂ @ 2.5 ppm v/v Acetone/AA if NO INF Correction




SIS Data: VOX in CO₂ @ 0.27 ppm v/v AA + Acetone INT Correction



SIS Data: VOX in CO₂ @ 0.1 ppm v/v AA + Acetone INF ON/OFF




SIS Data: VOX in CO₂ @ 0.07 ppm v/v AA – if No Acetone present




SIS Data: O₂ in CO₂ @ 16 ppm v/v (CO₂ MK-UP System Employed)

Use of CO₂-based Cal Gas Stds is recommended for all Analytes



SIS Data: H₂O in CO₂ @ 8 ppm v/v Level

Use of Indirect H₂O Cal Method @ SIG = Xe is viable option



SIS – Data for Other ISBT Listed CO₂ Impurities

| Impurity | SIG | AM FR = Frag IS = Isotope | LDR ppm v/v | DL ppm v/v | Interference (s) |
|------------------|-----|---------------------------------|----------------|---------------|--|
| O ₂ | Xe | 32 | 0 – 40+ | LT 1 | None found-to-date |
| H ₂ O | Xe | 18 | 0 – 30+ | LT 1 | None found to-date |
| NO | Hg | 30 | 0 - 5+ | LT 0.05 | None found-to-date |
| NO ₂ | Hg | 46 | 0 - 5+ | LT 0.05 | DME = INT AM = 46 |
| NH ₃ | Hg | 17 | 0 - 5+ | LT 0.05 | None found-to-date |
| HCN | Xe | 27 | 0 - 5+ | LT 0.1 | VCI = INT HCN Lower DL @ SIG = Xe + Kr |
| VCI | Hg | 64 ^{is} | 0 - 5+ | LT 0.05 | Use VCI Isotope AM = 64 |

Note 1: PH₃ study not completed SIG = Hg, AM = 34

Note 2: CO work was not successful for this ISBT application

SIS Applicability to ISBT CO₂ List

| <u>ISBT LISTED IMPURITY</u> | <u>SIS Measurable?</u> | <u>Comments</u> |
|--|--|--|
| CO₂ Purity | NO (but positive CO ₂ ID) | Can't do N ₂ , H ₂ , Ar |
| Moisture (H₂O) | YES (use CO ₂ Make-up Gas) | Can go sub-ppm range |
| Oxygen (O₂) | YES (use CO ₂ Make-up Gas) | Requires O ₂ in CO ₂ Stds |
| Carbon Monoxide (CO) | NO | Poor Sensitivity issues |
| Ammonia (NH₃) | YES | Needs passivated inlet |
| NO / NO₂ | YES (+ speciation) | Some minor INF's (not common) |
| Non-Volatile Residue (NVR) | NO | Can't do particulates / non volatiles |
| Non-Volatile Organic Residue (NVOR) | NO (but C6+ = Yes) | Can't do V. Hi MW oils/grease |
| Phosphine (PH₃) | YES (+ H ₂ S) PH ₃ = RARE | Can't distinguish from H ₂ S |
| Total Volatile Hydrocarbons (THC) & TNMHC | THC = YES TNMHC = YES | Get CH ₄ + Algorithm Sum of all TNMHC's vs THA-FID signal |
| Acetaldehyde (AA) | YES (+ VOX speciation) | Needs corrections for INF's |
| Aromatic Hydrocarbons (AHC) | YES (B-T-EX) | V good precision / low DL's no known INF's |
| Total Sulfur Content (TSC*) less SO₂ | YES (+ speciation) | Algorithm Sum of VSC's |
| Sulfur Dioxide (SO₂) | YES | Easily Speciated from COS |
| HCN | YES | V. Low DL no common INF |
| Vinyl Chloride | YES | V Low DL |

SUMMARY

SIS Performance vs Traditional Bev-Gas Analyzers

- **Speed**: Advantage vs many analyzer systems (e.g. LT 1 min vs 10-20 min for full ISBT analyte list)
- **Footprint**: Good vs common analyzer rack systems
- **Versatility**: Flexible list of monitored impurities vs most traditional Analyzers
- **Speciation Capability**: Less than GC but advantages vs single analyte analyzers. Minimizes DT use.
- **Measurement Quality**: LDR / Precision / DL: equivalent or superior to traditional analyzers
- **Freedom from Interferences**: Mixed – Application dependent: For AHC - is superior to many systems. Awareness of interferences is key to success (e.g. AA). Proper SIG / AM selection, Sequencing + Software-correction routines are effective in most cases.
- **Calibration**: Indirect Calibration = advantage. Recommend CO₂ based Cal gas – a disadvantage
- **Support**: Does not require He, N₂, H₂, Zero Air or expert/skilled operators
- **Maintenance**: Need for periodic SIG gas replacement + SIG Filament changes



Where do we go from Here?


R&D – New CO₂ SIS Methods Development Feedgas + in-Process Gas (e.g. Before/After Filter Beds = instant bed breakthrough warning) + Captured Recycle Gas use!




-Trace O₂ in Bev-Grade N₂ (NEW ISBT N₂ Guidelines coming)

- Impurities in PET bottles, preforms & resins (e.g. AA levels)


SIS-based Methods for trace aromatics in new / aged packaged beverages (shelf-life / stability testing)




-AND Beyond??




CO₂ Feedgas Method R&D - NatWells



CO₂ Feedgas Method R&D – NatWells Method Optimization In-Progress



CO₂ Feedgas Method R&D - BioFerm



Acknowledgements

Our Special Thanks to:




Mr. Rodney Resch & Mr. Walter
Woodland – *V&F* Instruments

Questions??



Airborne Labs International

SIS Analyzer – Under-the-Hood



Hg Reservoir

Xe Solenoid

SIG Dual Filament Ass'y


Hg Htr

Xe Inlet Line

Sample / Span Gas / Zero

Gas Inlet Manifold & Solenoids

Left Side View →




Turbo-Pump

Xe & Kr Cylinders

Hg C-Trap


The "Soft Ionization" Process — NEW CO₂ System



New Inlet Pressure Make-Up System Developed

3 Optimized ISBT Cal Gas Stds in CO₂ Developed

New CO₂ Make-up System



Required for accurate trace O₂ + H₂O SIS measurements

NEW - ISBT List SIS Calibration Strategy **(3 Combo Cal Std Cylinders – CO₂ Balance Gas + Certified – 3-stage filtered CO₂ Zero Gas Std)**

SIS Cal Std #1: Benzene – COS – AA – CH₄

SIS Cal Std #2: Acetone – DME – DMS – O₂

SIS Cal Std #3 Option 1: CO – Ethane – MeOH – NO₂
or

SIS Cal Std #3 Option 2: CO – n-Pentane – EtOH – NO₂

Certified Zero Gas – ISBT-List Tested 20 lb Tank of High Purity CO₂ + 3-bank Filtered (H₂O – Charcoal – O₂ scrubber cartridges)

List of key Cal Std analytes selected -- based on SIG Gas + AM's used + interference correction routines used & ability to act as accurate surrogate stds for other ISBT List target analytes (e.g. H₂S, SO₂, NO, TEX, H₂O, HCN, etc)



NEW THC / TNMHC / TCS*+ TOVO Algorithms

Required for SIS data correlation to existing (historical) THC (FID), GC and TSC Analyzer Data + ISBT / Bottler Impurity Guideline Limits

THC (FID) Equivalence: R_{ij} 's of all common THC (FID) active impurities (e.g. C2-C5 VHC's + VOX + VSC targets established & algorithm written – to convert SIS data into THC(FID) “equivalence as CH₄” values + TNMHC calculated (THC – CH₄).

TSC* Equivalence: SIS Algorithms written to properly calculate TSC* value (based on ISBT Method 14.0) of all SIS speciated VSC targets – less SO₂.

TOVO (Total “Other” Volatile Oxygenates) Equivalence: SIS Algorithm written to calc TOVO from SIS speciated VOX data.