

Passivated Gas Analysis Hardware

Extremely Favorable Pricing & Rapid, Worldwide Delivery

KEY BENEFITS

- Airborne Labs International (ALI) offers a diverse range of high quality passivated & stainless steel gas handling hardware at **extremely competitive** pricing! We also provide you with rapid **worldwide delivery** through our regional distributors. **Just compare to be aware!**
- Our hardware offerings are ideal for high purity liquefied gas sampling, gas analyzer systems, and challenging ppm/ppb level gas standards & samples involving “sticky” sulfur agents, or other challenging trace impurities.
- Our **Silconert® passivated** & stainless steel hardware products include: 1 & 2 Stage Mini-Pressure Regulators, Low Pressure Warning Sensors (**LoGas™**), CGA adapters, all types of valves (ex. ball, ball with vent, needle, needle with pressure reliefs, fine metering, toggle, check valves), long coils of transfer line tubing, flexible tubing, unions, tees, crosses, elbows, reducers, family adapters, special fixed purge flow assemblies, high pressure quick connects (with protective lanyard caps & low temperature compatible gaskets), in-line filters - **and more!**

PRODUCT BRIEF
PASSIVATED HIGH PURITY GAS ANALYSIS HARDWARE



1-Stage Passivated & SS Pressure Regulators



2-Stage Passivated & SS Mini-Pressure Regulators




LoGas™ Low Pressure Warning Sensor



Passivated & SS CGA Adapter Assemblies



Passivated & SS High Flow / Low Temp Quick Connects



Passivated & SS Transfer Line Tubing Coils



Passivated, High Pressure Flexible Tubing



Passivated & SS Metering Valves



Passivated & SS Needle Valves



Passivated & SS Needle Valves with Burst Disk



Passivated & SS Toggle Valves



Passivated Needle Valve - Burst Disk - Dip Tube Assemblies



Passivated & SS Fine Metering Valves



Passivated Fixed Bypass CV Assembly 3.2 LPM



Passivated & SS Ball Valves



Passivated & SS Venting Ball Valves



Passivated & SS 90 Elbows



Passivated & SS Check Valves



Passivated & SS Adapters



Passivated & SS Crosses



Passivated & SS Tee Fittings



Passivated & SS Unions

What do you need? - Yes we have that!



Airborne Labs International

ISO / IEC 17025:2005 Accredited Laboratory

AIRBORNE LABS INTERNATIONAL, INC
22C World's Fair Drive, Somerset, New Jersey 08873
T: (+1) 732-302-1950 F: (+1) 732-302-3035
E: Sales@Airbornelabs.com www.AirborneLabs.com