

# Sensory 15/16 Kit™

FOR PERFORMING ISBT METHODS 15 & 16 SENSORY TESTS

## KEY BENEFITS

- Performs ISBT Method 15 “Odor & Appearance of Snow” and ISBT Method 16 “Odor, Taste & Appearance in Water.”
- Conical flask and donut weight design ensure sampling and analysis safety. Special inlet & outlet cap assemblies allow for easy connection to a pressure limited sample source (ex. ZCheck-XLT).

## SENSORY KIT COMPONENTS

- 2x heavy wall 250 cc conical flasks
- 2x impinger heads
- 2 inlet / outlet special cap assemblies
- 2x donut weights
- 200x 3 ounce taste test cups
- ISBT Test Instructions

## RECOMMENDED ACCESSORIES

- **VC-Versipad / 1 Liter Sample Cylinder / Flash Vaporizer-Pressure Regulator (VPR)**  
Used for proper flash vaporization of a Liquid CO<sub>2</sub> sample and proper pressure delivery to a downstream flow controller (ex. ZCheck-XLT).
- **ZCheck-XLT**  
Highly recommended safety equipment used to limit pressure and accurately control sample flow to Sensory Kit 15/16 glassware. Includes a 5 psig check valve, DT-5 digital timer, passivated fine metering valve / toggle valve, beverage grade plastic tubing & special sleeve connector.

## RECOMMENDED ACCESSORIES



SENSORY 15/16™ KIT



SNOW POWDER FUNNEL



SNOW POWDER SPATULA



VC-VERSIPAD™ WITH 1 LITER SAMPLE CYLINDER AND VPR



SNOWQUIK 1.3B™



1L POLISHED NVR CAN



ST-6 SAMPLE TRANSFER LINE



FIBERLESS SYNTHETIC SNOW BAG



ZCHECK-5XLT™



**Airborne Labs International**

ISO / IEC 17025:2005 Accredited Laboratory

**AIRBORNE LABS INTERNATIONAL, INC**  
22C World's Fair Drive, Somerset, New Jersey 08873  
T: (+1) 732-302-1950 F: (+1) 732-302-3035  
E: Sales@AirborneLabs.com www.AirborneLabs.com