

True Blue™ MLB Gas Sampling Bags

ISBT-GRADE, CERTIFIED CLEAN GAS SAMPLING BAGS

KEY BENEFITS

- These highly rugged **True Blue™ Multi-Layer Barrier (MLB)** bags exhibit **excellent cleanliness** and **sample stability** properties. They meet or exceed all guidelines for sampling as described in **ISBT Method 1.2** and are superior to Tedlar® film bags in many applications.
- ALL's proprietary MLB film, valve materials and manufacturing process **eliminate gas permeability and bag material outgassing issues**.
- Each lot of MLB bags is analyzed for cleanliness and a **Certificate-of-Cleanliness** issued. Available in 2L or Custom Order bag sizes range from 1L to 20L.
- **MLB** bags are sold either individually or in cases of 10 bags.
- **MLB** bags are suitable to sample, store and ship as **non-hazmat** for all Division 2.2 gas samples including: CO₂, N₂, Ar, Air, SF₆ and many halocarbon-based refrigerants. **MLB** bags do not need hermetic air-shipping cans and can be filled to between 60-80+% capacity without leakage issues. Storage stability times for sample impurities typically **exceed 15 days**.
- Many important industrial, environmental and medical gas sampling applications can employ **MLB** bags to obtain accurate test results.

TYPICAL APPLICATIONS

- Gas Production Facilities
- Distribution Depot Facilities
- Bottling Plants
- Brewery Gas Testing
- Restaurant Fountain Minibulk Tanks
- QC Laboratory
- Power Plants



MLB BAG - 2L



SHA-316 T-FITTING



Airborne Labs International

ISO / IEC 17025:2005 Accredited Laboratory

AIRBORNE LABS INTERNATIONAL, INC

22C World's Fair Drive, Somerset, New Jersey 08873

T: (+1) 732-302-1950 F: (+1) 732-302-3035

E: Sales@AirborneLabs.com www.AirborneLabs.com