

SnowQuik™ 1.3 Series

LIQUID CO₂ SNOW GENERATORS

KEY BENEFITS

- Simple, Compact, Versatile, Rugged, Economical snowmaker that generates approximately 300+ g of snow in 1 minute for ISBT / CGA specified tests.
- A & C-Model is specifically designed for hands-free use in CGA or ISBT methods for non-volatile residue (NVR), non-volatile organic residue (NVOR) & sensory tests from Liquid CO₂.
- 6 ft. (3 m) x 1/4 inch stainless steel braided, flexible PTFE lined transfer line with 1900 psig (13,100 kPag) pressure relief device for operational convenience and safety. (Model B & C)
- Full flow high pressure / low temperature quick connect fittings allow for **tool-free** source connections.

MODELS

- **SnowQuik-1.3A**
Hard Mounted snow generator with quick connect cone for easy, convenient storage tank, tank truck, rail car, or minibulk tank sampling.
- **SnowQuik-1.3B**
Hand held snow generator with quick connect pressure relief assembly for easy, convenient storage tank, tank truck, rail car, or minibulk tank sampling.
- **SnowQuik-1.3C**
Trigger Clamp with quick connect pressure relief assembly inlet fitting for safe *Hands-Free* snow sample generation from storage tanks, tank trucks, or minibulk tanks.

TYPICAL APPLICATIONS

- ISBT Method 1.3
Carbon Dioxide Sampling Using Snow Generators and Containers
- ISBT Method 8.0
Total Non-Volatile and Non-Volatiles Organic Residues by Gravimetry
- ISBT Method 15.0
Odor of Carbon Dioxide Snow

SnowQuik-1.3A



SnowQuik-1.3B



SnowQuik-1.3C



RECOMMENDED ACCESSORIES



Digital Scale (DS-2KG)



Female QC w/ dust plug



1L Polished NVR Can



Cryo-Safety Gloves



Fiberless, Synthetic Snowbag
(Snowbag White)



ISO / IEC 17025:2017 Accredited Laboratory

©2010 Airborne Labs International, Inc. All rights reserved. The Airborne Labs International logo and design are registered trademarks of Airborne Labs International, Inc. All other trademarks not owned by Airborne Labs International, Inc. or its subsidiaries that are depicted herein are the property of their respective owners. Airborne reserves the right to change this document at any time without notice and disclaims liability for editorial, pictorial or typographical errors.

AIRBORNE LABS INTERNATIONAL, INC
22C World's Fair Drive, Somerset, New Jersey 08873
T: (+1) 732-302-1950 F: (+1) 732-302-3035
E: Sales@AirborneLabs.com www.AirborneLabs.com